



LS486 Tek-MMARS®400 Software IL790 Specification Sheet

Tek-CARE®400 and Tek-CARE®400 P5 Nurse Call Systems • Rev. 2 - 11/2012

Features

- Centralized and network-accessible call monitoring.*
- Centralized and network-accessible reporting via web browser.

* Not intended as the primary means of annunciation.

Description

TekTone®'s LS486 Tek-MMARS®400 Management Monitoring And Reporting Software provides real-time nurse/patient communications monitoring. Tek-MMARS®400 allows nurse call system activity to be monitored on the Tek-BRIDGE™ PC and stored for follow-up review and evaluation.

Reports are viewed using a web browser (such as Microsoft® Internet Explorer). They may also be printed, saved as HTML files for later viewing in a web browser, or saved as CSV files for viewing in a spreadsheet program (such as Microsoft® Excel).

When connected to a facility's LAN, the Tek-BRIDGE™ PC provides these services on networked computers via a web browser.


Listings



UL® 1069 listed and cUL® listed to CSA C22.2 No. 205; UL® File Number E71599.


Required Components

Tek-CARE®400 or Tek-CARE®400 P5 nurse call system
NC470 Tek-BRIDGE™ PC
LS486 Tek-MMARS®400 Software Component License
Optional connection to the facility's network



Event Monitor

Reporting



LS486 Tek-MMARS®400 Software

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)